

Product datasheet for **MR210216L3V**

Lrp10 (NM_022993) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lrp10 (NM_022993) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lrp10
Synonyms:	Lrp9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022993
ORF Size:	2142 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210216).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022993.3 , NP_075369.2
RefSeq Size:	3013 bp
RefSeq ORF:	2142 bp
Locus ID:	65107
UniProt ID:	Q7TQH7
Cytogenetics:	14 C2
Gene Summary:	Probable receptor, which is involved in the internalization of lipophilic molecules and/or signal transduction. May be involved in the uptake of lipoprotein APOE in liver. [UniProtKB/Swiss-Prot Function]



[View online »](#)