

Product datasheet for **MR210224L3V**

Mtss1l (NM_198625) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Mtss1l
Synonyms:	ABBA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_198625
ORF Size:	2145 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR210224).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_198625.1 , NP_941027.1
RefSeq Size:	4760 bp
RefSeq ORF:	2148 bp
Locus ID:	244654
UniProt ID:	Q6P9S0
Cytogenetics:	8 E1



Gene Summary:

Involved in plasma membrane dynamics. Potentiated PDGF-mediated formation of membrane ruffles and lamellipodia in fibroblasts, acting via RAC1 activation. May function in actin bundling.[UniProtKB/Swiss-Prot Function]