

Product datasheet for **MR210232**

Dpy19l3 (NM_178704) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Dpy19l3
Synonyms:	6030410G08; 9330164H19Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >MR210232 representing NM_178704
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGATGTACATCCGGCAAAGAAAAGAGACAAAACCCATAGAAGTTTCTGAGGATTTTCCATACCAAAAAG
 AAGATGTGAAGTTGAAAAGAAATTGCCATCTGGTTGTGCCTCTGGAAGTTCTGGAAGATCCTGTCATC
 CGCAGTTGGCGGAAGTGTGGCTCTTTGCATTGGACTTCTTACATCCGTTTACCTGGCCACCTGCATGAG
 AACGATTTGTGGTTTTCCAATATTAAGGAAGTAGAGCGAGAGATCTCATTACAGGACGGAATGTGGACTGT
 ACTATTTCTTACTACAAGCAGATGCTGCAGGCTCCGACCCTCCTGCAAGGTTTCCATGGCTTAATCTATGA
 CAATAAAAGTGAATCCATGAGGACGATCAACCTCCTTCAACGGATGAACATTTACCAAGAAGTTTTTCTC
 AGTGTTTTATATAGAGTTCTCCCCATACAGAAATATCTAGAACCTGTTTATTTTATATCTACACTCTGT
 TCGGGCTCCAAGCAGTCTATGTCACGGCCCTTATATAACCAGCTGGCTCCTAGTGGCACCTGGCTGTC
 GGGACTGTTGGCAGCCCTCTGGTACGTACCAACAGAAATAGACACCACAAGAGTCGAGTTTACCATCCCG
 CTGAGGGAGAACTGGGCTCTGCCTTTTTTGCATTCAGATTGCAGCAATTACCTACTTCTGAGACCAA
 ATTTACAGCCTCTTTCTGAAAGGCTGACACTTCTTGCCATTTTCGTATCAACTTTTCTTTAGCCTGAC
 ATGGCAATTTAACCAGTTTATGATGCTGCTGCAGGCTGGTGCTTTCATCCTGGACTCCCTGGACATG
 CTGCCGGCCATGAAGGCCACCTGGCTGTACGGGATACAGATAAGCTGCCTGCTCCTCGTCTGCACGCTTC
 AGTTCTTTAACTCCATGATTCTGGGATCGTTGCTCATCAGTTTTAACCTTTCAGTATTAATTGTAAGAAA
 ACTTCAGAAAAACCTGAAGACTGGATCCTTCTGACTAGAATTTGAAACTACTGTTGCATTTGCTTCTA
 GTTTTCTGTTGACACTTTTTCTCAACAACATTATTAAGAAAGTTCTTAACCTGAAGTCAGATGAACACA
 TCTTTAAATTTCTGAAGCCAAATTTGGGTTTGGAGCAACAAGGGATTTTCGATGCAATCTCTATCTGTG
 TGAGGAGGCCCTTCGGCCTCCTGCCTCTCAACACGTTCCAGAGGCTCTCGGAGACACTGCTTTTTTACGCA
 TACATGTTTGTCTCGTGGTCACAGTGGTTACAGCATCCGTTGTAGCCTTTTATAATCTCAGTGATTCCA
 CAAGTCTGAAGTCCATGGATCAGACAAGAAAGCGTGCTGTTGACCTGAAGCCAGAAGCTGCGTACAACCT
 GATCCATACCATTTCTGTTGGGGTCTTGGCGTTGAGTACGATGAGGATGAAGTACCTGTGGACCTCCAC
 ATGTGTGTGTTTCGGCTCCTTCGGCCTGTGCAGCTCTGAAGTCTGGGAAGTGTCTACTGAGGCTGGTTTCATC
 TCTGTAACCCAAAGCGGATCTGGGTCTGCGCTACCTGGTGCCGTTCTCACCTGCTGTACCTGTGCTA
 TAAGTCTGGCCAGGAGTATGGATGAAGTCTCGGAGCTGAAAGAATTCTACGACCCGGATACAGTGGAG
 CTGATGACCTGGATTAACTCTAACACACCAAGGAAGGCTGTGTTGCTGGAAGCATGCAGTTGCTGGCTG
 GAGTCAAGCTGTGCACAGGACGACCTAACCAACCCACACTATGAAGATAAAAGTCTGAGAGAGCG
 GACCCAGGCGGTTATCAGATCTATGCGAAGAGGTGCGCAGAGGAAGTGCACGCACTGCTGAGGTCCTTT
 GGCACCTGACTTCGTGATCCTGGAGGACAGCATCTGCTACGAACGACGACACAGGGGCTGCCGGCTGC
 GCGACTTGCTGGATGTGGCCAACGGACATGAAATGGATGGGCCAGGAGAGAGTGACCTGACTTGAGACC
 TGCAGACCAACCTCGCTTCTGTGAAGAGATTAAGGAACCTGCCCTCCTATGCAGCCCACTTCACTAGA
 GTATTCAGAAACAAGACATTCCATGTCTACAAGCTGTCCAGAAACAAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR210232 representing NM_178704

Red=Cloning site Green=Tags(s)

MMYIRQRKETKPIEVSEDFPSKEDVKLEKKLPSGCASGRFWKILSSAVGGTVALCIGLLTSVYLATLHE
NDLWFSNIKEVEREISFRTECLYYSYKQMLQAPTLQGFHGLIYDNKTESMRTINLLQRMNIYQEVFL
SVLYRVLPIQKYLEPVYFYIYTLFGLQAVYVTALYITSWLLSGTWL SGLLAALWYVTNRIDTTRVEFTIP
LRENWALPFFAIQIAAITYFLRPNLQPLSERLTLLAIFVSTFLSLTWQFNQFMMLLQALVLFILDSLDM
LPAMKATWLYGIQISCLLLVCTLQFFNSMILGSLLSIFNLSVLIVRKLQKNLKTGSFLTRIWKLLHLL
VFCLTLFLNNIIKKVLNLKSDEHIFKFLKAKFGFGATRDFDANLYLCEEAFGLPLNTFQRLSETLLFYA
YMFVLVVTVTASVVAFHNLSDSTSLKSMQTRKRAVDLKPEAAYNL IHTILFGVALSTMRM KYLWTS
MCVFASFGLCSSEVWELLLRLVHLCNPKRIWVLRYLVPVLTLLYCYKSWPGVMDELSELKEFYDPDTVE
LMTWINSNTPRKAVFAGSMQLLAGVKLCTGRTL TNHPHYEDKSLRERTQAVYQIYAKRSPEEVHALLRSF
GTDFVILEDSICYERRHHRGCLRDLDDVANGHEMDGPGESDPDLRPADHPRFCEEIKRNLPSYAAHFTR
VFQNKTFHVYKLSRNK

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/mm9027_c01.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:


NM_178704

ORF Size:

2148 bp

OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery. The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_178704.3 , NP_848819.2
RefSeq Size:	7671 bp
RefSeq ORF:	2151 bp
Locus ID:	233115
UniProt ID:	Q71B07
Cytogenetics:	7 B2
MW:	83.4 kDa
Gene Summary:	Probable C-mannosyltransferase that mediates C-mannosylation of tryptophan residues on target proteins.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR210232