

Product datasheet for **MR210249L2V**

Eif2b5 (NM_172265) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Eif2b5 (NM_172265) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Eif2b5
Synonyms:	C81315
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_172265
ORF Size:	2151 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210249).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172265.2 , NP_758469.1
RefSeq Size:	2556 bp
RefSeq ORF:	2154 bp
Locus ID:	224045
UniProt ID:	Q8CHW4
Cytogenetics:	16 A3
Gene Summary:	Catalyzes the exchange of eukaryotic initiation factor 2-bound GDP for GTP.[UniProtKB/Swiss-Prot Function]



[View online »](#)