

Product datasheet for **MR210256L4V**

Sfil (BC046305) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|----------------------------------|-----------------------------------|
| Product Type: | Lentiviral Particles |
| Symbol: | Sfil |
| Synonyms: | mKIAA0542, RP23-191E3.10 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | BC046305 |
| ORF Size: | 2151 bp |

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR210256).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [BC046305](#), [AAH46305](#)

RefSeq Size: 3722 bp

RefSeq ORF: 2153 bp

Locus ID: 78887

Cytogenetics: 11 A1

Gene Summary: Plays a role in the dynamic structure of centrosome-associated contractile fibers via its interaction with CETN2.[UniProtKB/Swiss-Prot Function]



