

Product datasheet for **MR210266L4V**

Lrfn5 (BC052038) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lrfn5 (BC052038) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lrfn5
Synonyms:	AI427653; AI604817; C130061B21; mKIAA4208
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC052038
ORF Size:	2157 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210266).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC052038 , AAH52038
RefSeq Size:	2782 bp
RefSeq ORF:	2159 bp
Locus ID:	238205
Cytogenetics:	12 C1
Gene Summary:	Cell adhesion molecule that mediates homophilic cell-cell adhesion in a Ca(2+)-independent manner. Promotes neurite outgrowth in hippocampal neurons (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)