

Product datasheet for **MR210347L4V**

Lss (NM_146006) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lss (NM_146006) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lss
Synonyms:	2810025N20Rik; D10Ertd116; D10Ertd116e; O; Osc
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_146006
ORF Size:	2199 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210347).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_146006.2 , NP_666118.1
RefSeq Size:	3192 bp
RefSeq ORF:	2202 bp
Locus ID:	16987
UniProt ID:	Q8BLN5
Cytogenetics:	10 39.1 cM



[View online »](#)

Gene Summary:

This gene encodes an enzyme that catalyzes the conversion of (S)-2,3 oxidosqualene to lanosterol. The encoded protein is a member of the terpene cyclase/mutase family and catalyzes the first step in the biosynthesis of cholesterol, steroid hormones, and vitamin D. [provided by RefSeq, Feb 2015]