

Product datasheet for **MR210353L3V**

Mthfd1l (BC030437) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mthfd1l (BC030437) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mthfd1l
Synonyms:	2410004L15Rik; AI647056; Fthfsdc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC030437
ORF Size:	2199 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210353).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC030437 , AAH30437
RefSeq Size:	3061 bp
RefSeq ORF:	2201 bp
Locus ID:	270685
Cytogenetics:	10 A1
Gene Summary:	May provide the missing metabolic reaction required to link the mitochondria and the cytoplasm in the mammalian model of one-carbon folate metabolism in embryonic an transformed cells complementing thus the enzymatic activities of MTHFD2.[UniProtKB/Swiss-Prot Function]



[View online »](#)