

Product datasheet for **MR210377**

Ap4b1 (NM_001163552) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ap4b1 (NM_001163552) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ap4b1
Synonyms:	1810038H16Rik; Ap4b4; AV004952
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>MR210377 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCATACCTTGGCTCTGAGGACGTGGTGAAGGAACTGAAGAAGGCTCTGTGTAACCTCATATTCAGG
 CTGATAGGCTGCGCTACCGAATGTCATCCAGCGAGTTATTAGGCACATGACTCAGGCTTGGACATGTC
 TGATGTTTTTCATGAAAATGGTAAAGGCCAGTGCCACTGTAGATATTGTTCAAAAAAGTTGGTTTATCTG
 TATATGGGCACATACGCACCTTTAAAACAGATCTAGCACTCTTGGCCATCAATACGCTCTGCAAAGACT
 GCTCAGACCCCAACCAATGGTGCAGGACTGGCACTACGGAGTATGTGTAGCCTCAGGATGCCTGGTGT
 GCAGGAATATATCCAGCAGCCTGTTGTTAACGGTCTGCGTGATAAAGCTTCATATGTCAGAAGAGTAGCA
 GTCCTCGGCTGTCCAAAATGCATAATCTTCACGGAGACTCTGAAGTGGATGGTGTCTCGTCAATGAGT
 TGTACAGTCTTCTGCGTGACCAGGATCCAATTGTGGTTGTGAAGTGTCTGAGGTCTCTAGAGGAAATTTT
 AAAGCAGGAAGGAGGTGTTGTCATTAATAAGCCCATTGCCACCATCTCTAAATAGAATGTCAAAGCTG
 GACCAGTGGGGCCAGGCTGAAGTACTGAACTTTCTGCTCCGCTACCAACCCCGCAGTGAGGAGGAACTGT
 TTGATATTCTCAATTTGTTGGACAGCTATCTCAAGAGCAGCAGTACTGGCGTGGTGTGAGGACCCAA
 ATTCTTTCTGATTTTGGCAAAAAGTTCACCCACGTGCAGACTGATGTGCTTGTGCGAGTCAAGGGGCC
 TTGCTGGCCGCTGCTCCTCAGAGAGCCGTGAGCTCTGTTTTGCTGCTCTCTGCCATGTCGCCCAGGTCT
 TGCACAGTTTGGCAGGCTCACTTTCAGCAGTCACTACAAAAGTTTTTTTGTCTTACTCGGAACCCCACTA
 CATCAAGCTCCAGAAGGTGGAGGTGCTGTGCGAGCTGGTGAATGACGAGAACGTGCAGCAGGTGCTGGAA
 GAGTGCCTGGTACTGCACTGACGTGGCTGCTGACTTTGCACAGGCTGCCATCTTCGCCATAGGTAGCA
 TTGCCAAGACTTACACAGATCAATGTGTGCAGATTCTAACAGAGTTGCTGGGACTTCGACAAGAACATAT
 CACTACTGTGGTTCGTGCAGACTTTCCGAGACCTGGTCTGGCTATGCCCTCAGTGTACTGAAGCTGTGTGT
 CAGGCTCTGCCCGCTGTGAGGAGAGCATTCAAGACAGTGAAGGAAACAGGCACTTATTTGTTACTTG
 GGGTCCACGGGAAAAAATCCCCAATGCTCCTTATGTGTTGGAGGACTTTGTAGACAACGTGAAGTCAGA
 GACATTTCCCGCTGTTAAGATGGAGTACTGACCGCACTGATGCGCCTTGTCTCTCGCGCCCGCTGAG
 TGCCAGGACATGCTGGGACGGCTGCTGCATTACTGCATAGAGGAAGAAAAGGATATGGCTGTACGGGACC
 GGGGTCTCTTCTATTACCGCCTCTGTTAGTTGGCATTGATAAAGTTAAGCAGATCCTGTGTAGTCTTAA
 GTCTGACCCTTCTTTGGACTTTTGGAGGATCAGCCGAAACGACCAGTAAATAGCTGGGCTTCAGACTTC
 AACACACTGGCGCCAGTGTATGGCAGAGCCCACTGGGCAACTATCTCTAAATGCCAGCAGGTAGAGCGGC
 ATCGCCTTGAGCTTCCCCACAATGCATCCTTTGCCACATCAGGGCACTTGATTTCCGAGGAGAACAAAA
 AGGTGCACAAGAACCTGATTCTGCACTCTCATGCTGGTCCCTAATCTCCAGCTTACTGTGAGTATTTT
 GAGAAGACGTGGCTTAGCCTCAGAGTTAGTTACCAGCAGGTGTTTCTTGGCAGGGAGAAGTCCAGCCTG
 ACACTCTACAAATGGCTCTAAAAGTAGTGAACATTCAGACCATCGCAATGAGCAGAGCTGGGGCACAGCC
 CTGGAAGCCTACCTCAGTGTCTCAGGATGATACTGGCGGCCTTTTCTTAGCTGAATTGCTACTGAAGCCT
 GAGAACTCAGAGATGCAATCTCTGTGAAGCAAAGTAAGGCAAGGACTGAATCCCTGCATGGTTTTGTTT
 CTGTCTGGAACTGTGATTGGAACAGTTGGAGACATAAAATCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >MR210377 protein sequence
Red=Cloning site Green=Tags(s)

MPYLGSEDEVVKELKKALCNPHIQADRLRYRNVIRVIRHMTQGLDMSDFMEMVKASATVDIVQKKLVYL
 YMGTYAPLKPDLALLAINTLCKDCSDPNPMVRGLALRSMCSLRMPGVQEYIQPVVNGLRDKASYVRRVA
 VLGCAKMHNLHGDSEVDGALVNELYSLLRDQDPIVVVNLRSLEEILKQEGGVVINKPIAHLLNRMSKL
 DQWGQAEVLNFLRLRYQPRSEEEELFDILNLLDSYLKSSSTGVVMGATKFFLILAKKFPHVQTDVLRVKG
 LLAACSSSERELCAALCHVRQVLHSLPGHFSSHYKKFFCSYSEPHYIKLQKVEVLCENVQVLE
 ELRGYCTDVAADFQAQAAIFAIGSIAKTYDQCQVQILTELLGLRQEHITTVVVQTFRDLVWLCPCQTEAVC
 QALPGCEESIQDSEKQALIWLLGVHGEKIPNAPYVLEDFVDNVKSETFPAVKMELLTALMRLVLSRPAE
 CQDMLGRLLHYCIEEEKDMAVRDRGLFYRLLLVGIDKVKQILCSPKSDPSLGLLEDQPERPVNSWASDF
 NTLAPVYGRAHWATISKCCQVERHRELPNASFATSGHLISEENKKAQEPDSDLMLVPLQLTAEYF
 EKTWLSLRVSYQQVFPWQGEVQPDTLQMALKVVNIQTIAMSRAGAQPWKAYLSAQDDTGGLFLAELLKPK
 ENSEMQISVKQSKARTESLHGFVSVLETVIGTVGDIKS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001163552

ORF Size: 2217 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001163552.1](#), [NP_001157024.1](#)

RefSeq Size: 2701 bp

RefSeq ORF: 2217 bp

Locus ID: 67489


UniProt ID: [Q9WV76](#)

Cytogenetics: 3 F2.2

MW: 83 kDa

Gene Summary: Component of the adaptor protein complex 4 (AP-4). Adaptor protein complexes are vesicle coat components involved both in vesicle formation and cargo selection. They control the vesicular transport of proteins in different trafficking pathways. AP-4 forms a non clathrin-associated coat on vesicles departing the trans-Golgi network (TGN) and may be involved in the targeting of proteins from the trans-Golgi network (TGN) to the endosomal-lysosomal system (By similarity). It is also involved in protein sorting to the basolateral membrane in epithelial cells and the proper asymmetric localization of somatodendritic proteins in neurons (PubMed:18341993). AP-4 is involved in the recognition and binding of tyrosine-based sorting signals found in the cytoplasmic part of cargos, but may also recognize other types of sorting signal (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR210377