

## Product datasheet for MR210400L3V

## Mfn1 (NM\_024200) Mouse Tagged ORF Clone Lentiviral Particle

## **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Mfn1 (NM_024200) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mfn1
Synonyms:	2310002F04Rik; 6330416C07Rik; D3Ertd265e; HR2; mKIAA4032
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024200
ORF Size:	2226 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210400).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<u>NM 024200.2, NP 077162.1</u>
RefSeq Size:	4531 bp
RefSeq ORF:	2226 bp
Locus ID:	67414
UniProt ID:	<u>Q811U4</u>
Cytogenetics:	3 15.75 cM



View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn



Gene Summary:Mitochondrial outer membrane GTPase that mediates mitochondrial clustering and fusion<br/>(PubMed:12527753, PubMed:23921378, PubMed:24513856, PubMed:15297672). Membrane<br/>clustering requires GTPase activity (By similarity). It may involve a major rearrangement of<br/>the coiled coil domains (PubMed:15297672). Mitochondria are highly dynamic organelles, and<br/>their morphology is determined by the equilibrium between mitochondrial fusion and fission<br/>events (PubMed:12527753). Overexpression induces the formation of mitochondrial<br/>networks (in vitro). Has low GTPase activity (By similarity).[UniProtKB/Swiss-Prot Function]

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US