

Product datasheet for **MR210425L3V**

Adarb2 (NM_052977) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adarb2 (NM_052977) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Adarb2
Synonyms:	Adar3; RED2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_052977
ORF Size:	2238 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210425).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_052977.3 , NP_443209.2
RefSeq Size:	5633 bp
RefSeq ORF:	2238 bp
Locus ID:	94191
UniProt ID:	Q9JL20
Cytogenetics:	13 A1
Gene Summary:	Lacks editing activity. It prevents the binding of other ADAR enzymes to targets in vitro, and decreases the efficiency of these enzymes. Capable of binding to dsRNA but also to ssRNA (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)