

Product datasheet for **MR210485L4V**

Exoc3 (NM_177333) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Exoc3 (NM_177333) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Exoc3
Synonyms:	2810050O03Rik; E430013E20Rik; Sec6; Sec6l1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_177333
ORF Size:	2268 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210485).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177333.3
RefSeq Size:	4597 bp
RefSeq ORF:	2268 bp
Locus ID:	211446
UniProt ID:	Q6KAR6
Cytogenetics:	13 C1
Gene Summary:	Component of the exocyst complex involved in the docking of exocytic vesicles with fusion sites on the plasma membrane.[UniProtKB/Swiss-Prot Function]



[View online »](#)