

Product datasheet for **MR210490**

Trpv2 (NM_011706) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Trpv2 (NM_011706) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Trpv2
Synonyms:	GRC; OTRPC2; VRL-1; Vrl1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR210490 representing NM_011706
 Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGACTTCAGCCTCCAACCCCCAGCTTTTAGGCTGGAGACGTCCGATGGAGATGAAGAAGGCAGTGCTG
 AGGTGAACAAAGGAAAGAATGAGCCACCCCATGGAGTCTCCCTTCCAGGGAGAGGACCGAACTTCTC
 CCCTCAGATTAAGTGAATCTCAACTACCGAAAGGGACTGGGTCCCAGCCAACAGGACCCAAATCGGTTT
 GACCGTGACCGACTCTTCAGTGTGGTCTCCCGGGGTGTCCCGAGGAGCTGACTGGACTGCTAGAGTACC
 TGCGCCGACCAGCAAGTACCTCACTGACTCGGCATACACAGAAGGCTCCACTGGAAGAGCTGCCTGAT
 GAAGGCTGTGCTGAACCTTCAGGATGGGGTCAATGCCTGTATCCTGCCGTGCTGCAGATTGACAGGGAT
 TCCGGCAATCCTCAGCCCTTGCAATGCCAGTGCACCGATGAGTTCTACCGAGGCCACAGTGCCTGC
 ACATCGCCATAGAGAAGAGGAGCCTGTGGTGCCTGAAGCTGCTGGTAGAGAATGGAGCGAATGTTACAT
 CCGAGCCTGTGGCCGCTTCTTCCAAAAGCACCAAGGAACCTGTTTCTATTTTGGAGAGCTACCTCTTTCT
 CTGGCAGCGTGCACCAAGCAGTGGGATGTGGTGCCTACCTCCTGGAGAACCACACCAGCTGCCAGCC
 TGGAGGCCACCGACTCCCTGGGCAACACAGTCTGCATGCTCTGGTAATGATTGCAGACAACCTCACTGA
 GAACAGTGCCTGGTGTATCCACATGTATGACAGCCTTCTCAAATGGGGGCCCGCTCTGCCCACTGTA
 CAGCTTAGGATATCTGCAACCATCAAGGCCTCACACCCCTGAAGTTGGCTGCCAAGGAAGTAAATTTG
 AGATCTTCAGGCACATCCTGCAGCGGGAGTTCTCAGGGCTGTACCAGCCCTTTCCCGAAAGTTCACCGA
 GTGGTGTACGGTCTGTCCGAGTGTACTGTACGACCTGTCTCTGTGGACAGTTGGGAAAAGAAGTCCG
 GCCTGGAAATCATCGCTTTCCATTGCAAGAGCCCGCACCGGCACCGCATGGTGGTTTTAGAGCCACTGA
 ACAAGCTTCTGCAGGAGAAATGGGATCGGCTCATCCCAAGATTCTTCTCAACTTCGCCTGTTACTTGGT
 CTACATGATCATCTTCACCATAGTTGCCTACCACAGCCTTCCCTGGAGCAGCCAGCCATTCCCTCATCA
 AAAGCGACTTTTGGGACTCCATGCTGCTGTTGGGCCACATTCTGATCCTGCTTGGGGGATTTTACCTCT
 TACTGGCCAGCTGTGGTACTTTTGGCGCGCGCCTGTTTCTGATCTGGATCTCGTTCATGGACAGTTACTT
 TGAAATCCTCTTCTTGTCCAGGCGCTGCTCACAGTGTGTCCCAGGTGCTGCGCTTCGTGGAGACTGAA
 TGGTACCTCCCCCTGTTAGTGTATCCCTAGTGTGGGCTGGCTGAACCTGCTTTATTATACACGTGGCT
 TTCAGCACACAGGCATCTACAGTGTATGATCCAAAAGGTCATTCTGCGAGACCTGCTCCGCTTCTGCT
 GGTCTACCTAGTCTTCTTTTCCGCTTTGCTGTAGCCCTAGTAAGCTTGAGCCGGGAGGCCCGAAGTCCC
 AAAGCCCCTGAAGATAGCAACCACTGTGACGAAAAGCCACGCTGGGTGAGGAGGAGCCAGTCC
 CATATGGGGGCATTCTGGATGCCTCCCTAGAGCTGTTCAAGTTCACCATGGTATGGGTGAGCTGGCTTT
 CCAGGAACAGCTGCGCTTTCGTGGGGTGGTGTGCTGTTGCTGTTGGCCTACGTCCTCCTCACCTACGTC
 CTAAGTGTCAACATGCTCATTGCCCTCATGAGTGTGAGACTGTCAACAGCGTTGCCACTGACAGCTGGAGCA
 TCTGGAAGTTGCAGAAAGCCATCTGTCTTGGAGATGGAGAAATGGTTACTGGTGGTGCAGGAGGAAAAG
 GCATCGCGCAGGGAGGCTGCTGAAAGTTGGCACCAAGGGGATGGTATACCTGATGAGCGCTGGTGTCTC
 AGGGTGGAGGAAGTAAACTGGGCTGCATGGGAGAAGACCCCTTCCACCTTATCTGAGGATCCATCAGGGG
 CAGGCATCACTGGTTATAAAAAGAACCCAACTCTAAACCTGGGAAGAACAGTGCCTCAGAGGAAGACCA
 TCTGCCTCTTTCAGTCTCCAGTCCCAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	<p>This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.</p>
Components:	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<p>NM_011706.2, NP_035836.2</p>
RefSeq Size:	<p>2751 bp</p>
RefSeq ORF:	<p>2271 bp</p>
Locus ID:	<p>22368</p>
UniProt ID:	<p>Q9WTR1</p>
Cytogenetics:	<p>11 B2</p>
MW:	<p>86.4 kDa</p>
Gene Summary:	<p>Calcium-permeable, non-selective cation channel with an outward rectification. Seems to be regulated, at least in part, by IGF-I, PDGF and neuropeptide head activator. May transduce physical stimuli in mast cells. Activated by temperatures higher than 52 degrees Celsius; is not activated by vanilloids and acidic pH.[UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for MR210490