

## Product datasheet for **MR210495L3V**

### Dact2 (NM\_172826) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Dact2 (NM_172826) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dact2
Synonyms:	2900084M21Rik; A630024E20; dapper2; Dpr2; Frd2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_172826
ORF Size:	2271 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210495).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_172826.3</a> , <a href="#">NP_766414.3</a>
RefSeq Size:	2849 bp
RefSeq ORF:	2274 bp
Locus ID:	240025
UniProt ID:	<a href="#">Q7TN08</a>
Cytogenetics:	17 A2



[View online »](#)

**Gene Summary:**

Involved in regulation of intracellular signaling pathways during development. Negatively regulates the Nodal signaling pathway, possibly by promoting the lysosomal degradation of Nodal receptors, such as TGFBR1. May be involved in control of the morphogenetic behavior of kidney ureteric bud cells by keeping cells epithelial and restraining their mesenchymal character. May play an inhibitory role in the re-epithelialization of skin wounds by attenuating TGF-beta signaling.[UniProtKB/Swiss-Prot Function]