

## Product datasheet for **MR210509L3V**

### Shroom1 (BC005761) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Shroom1 (BC005761) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Shroom1
Synonyms:	RP23-188H3.7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC005761
ORF Size:	2280 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210509).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC005761.1</a>
RefSeq Size:	2663 bp
RefSeq ORF:	2282 bp
Locus ID:	71774
Cytogenetics:	11 B1.3
Gene Summary:	May be involved in the assembly of microtubule arrays during cell elongation. [UniProtKB/Swiss-Prot Function]


[View online »](#)