

Product datasheet for **MR210523L4V**

Cep85 (NM_144527) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cep85 (NM_144527) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cep85
Synonyms:	2410030J07Rik; AI173272; AI173504; Ccdc21
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_144527
ORF Size:	2283 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210523).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144527.3 , NP_653110.3
RefSeq Size:	3878 bp
RefSeq ORF:	2286 bp
Locus ID:	70012
UniProt ID:	Q8BMK0
Cytogenetics:	4 D3
Gene Summary:	Acts as a negative regulator of NEK2 to maintain the centrosome integrity in interphase. Suppresses centrosome disjunction by inhibiting NEK2 kinase activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)