

Product datasheet for **MR210535L3V**

Cpxm2 (NM_018867) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cpxm2 (NM_018867) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cpxm2
Synonyms:	4632435C11Rik; CPX-2; Cpx2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018867
ORF Size:	2295 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210535).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018867.3
RefSeq Size:	3506 bp
RefSeq ORF:	2295 bp
Locus ID:	55987
UniProt ID:	Q9D2L5
Cytogenetics:	7 F3
Gene Summary:	May be involved in cell-cell interactions.[UniProtKB/Swiss-Prot Function]



[View online »](#)