

Product datasheet for **MR210590L4V**

SpeccII (BC037464) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	SpeccII
Synonyms:	mKIAA0376, cytospin A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC037464
ORF Size:	2328 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR210590).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	BC037464.1
RefSeq Size:	5181 bp
RefSeq ORF:	2330 bp
Locus ID:	74392
Cytogenetics:	10 C1



Gene Summary:

Involved in cytokinesis and spindle organization. May play a role in actin cytoskeleton organization and microtubule stabilization and hence required for proper cell adhesion and migration (By similarity).[UniProtKB/Swiss-Prot Function]