

Product datasheet for **MR210616L4V**

Aco2 (NM_080633) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Aco2 (NM_080633) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Aco2
Synonyms:	Aco-2; Aco3; D10Wsu183e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_080633
ORF Size:	2340 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210616).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080633.2 , NP_542364.1
RefSeq Size:	2785 bp
RefSeq ORF:	2343 bp
Locus ID:	11429
UniProt ID:	Q99K10
Cytogenetics:	15 E1
Gene Summary:	Catalyzes the isomerization of citrate to isocitrate via cis-aconitate.[UniProtKB/Swiss-Prot Function]



[View online »](#)