

Product datasheet for MR210640L4

Usp1 (NM_146144) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Usp1 (NM_146144) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Usp1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210640).
Restriction Sites:	Sgfl-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_146144
ORF Size:	2352 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_146144.3
RefSeq Size:	3567 bp
RefSeq ORF:	2355 bp
Locus ID:	230484
UniProt ID:	Q8BJQ2
Cytogenetics:	4 C6
Gene Summary:	This gene encodes a member of the ubiquitin-specific peptidase family. The encoded protein acts as a catalytic subunit in a heterodimeric deubiquitinating enzyme complex that deubiquitinates Fanconi anemia, complementation group D2, and plays a role in homologous recombination-mediated DNA repair. Disruption of this gene is associated with a Fanconi anemia-like phenotype and genomic instability. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene have been defined on chromosomes 3, 12, and 15. [provided by RefSeq, Aug 2014]

