

Product datasheet for **MR210646L3V**

Lrrc41 (NM_153521) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lrrc41 (NM_153521) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lrrc41
Synonyms:	AA409966; AW555107; D630045E04Rik; D730026A16Rik; MUF1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153521
ORF Size:	2358 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210646).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153521.1 , NP_705741.2
RefSeq Size:	2883 bp
RefSeq ORF:	2424 bp
Locus ID:	230654
UniProt ID:	Q8K1C9
Cytogenetics:	4 D1
Gene Summary:	Probable substrate recognition component of an ECS (Elongin BC-CUL2/5-SOCS-box protein) E3 ubiquitin ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.[UniProtKB/Swiss-Prot Function]


[View online »](#)