

Product datasheet for **MR210654L4V**

Itgb6 (NM_001159564) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Itgb6 (NM_001159564) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Itgb6
Synonyms:	2210409C20Rik; 4831415H04Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001159564
ORF Size:	2361 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210654).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001159564.1 , NP_001153036.1
RefSeq Size:	4894 bp
RefSeq ORF:	2364 bp
Locus ID:	16420
UniProt ID:	Q9Z0T9
Cytogenetics:	2 34.81 cM



[View online »](#)

Gene Summary:

Integrin alpha-V:beta-6 (ITGAV:ITGB6) is a receptor for fibronectin and cytotactin (By similarity). It recognizes the sequence R-G-D in its ligands (PubMed:10025398). ITGAV:ITGB6 acts as a receptor for fibrillin-1 (FBN1) and mediates R-G-D-dependent cell adhesion to FBN1 (By similarity). Integrin alpha-V:beta-6 (ITGAV:ITGB6) mediates R-G-D-dependent release of transforming growth factor beta-1 (TGF-beta-1) from regulatory Latency-associated peptide (LAP), thereby playing a key role in TGF-beta-1 activation (PubMed:10025398). [UniProtKB/Swiss-Prot Function]