

## Product datasheet for **MR210674L4V**

### Sf3a1 (NM\_026175) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sf3a1 (NM_026175) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sf3a1
Synonyms:	1200014H24Rik; 5930416L09Rik; A1159724
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026175
ORF Size:	2373 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210674).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026175.5</a> , <a href="#">NP_080451.4</a>
RefSeq Size:	4916 bp
RefSeq ORF:	2376 bp
Locus ID:	67465
UniProt ID:	<a href="#">Q8K4Z5</a>
Cytogenetics:	11 A1



[View online »](#)

**Gene Summary:**

Involved in pre-mRNA splicing as a component of the splicing factor SF3A complex that contributes to the assembly of the 17S U2 snRNP, and the subsequent assembly of the pre-spliceosome 'E' complex and the pre-catalytic spliceosome 'A' complex. Involved in pre-mRNA splicing as a component of pre-catalytic spliceosome 'B' complexes.[UniProtKB/Swiss-Prot Function]