

## Product datasheet for **MR210684L3V**

### Zc3h11a (BC066163) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Zc3h11a (BC066163) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zc3h11a
Synonyms:	mKIAA0663, Zc3hdc11a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC066163
ORF Size:	2376 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210684).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC066163.1</a>
RefSeq Size:	2815 bp
RefSeq ORF:	2378 bp
Locus ID:	70579
Cytogenetics:	1 E4
Gene Summary:	Involved in nuclear mRNA export; probably mediated by association with the TREX complex. [UniProtKB/Swiss-Prot Function]



[View online »](#)