

Product datasheet for **MR210694L3V**

Kifc3 (NM_001145832) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kifc3 (NM_001145832) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kifc3
Synonyms:	AI325457; BB123200
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145832
ORF Size:	2382 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210694).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001145832.1 , NP_001139304.1
RefSeq Size:	2945 bp
RefSeq ORF:	2064 bp
Locus ID:	16582
UniProt ID:	O35231
Cytogenetics:	8 47.12 cM
Gene Summary:	Minus-end microtubule-dependent motor protein. Involved in apically targeted transport. Required for zonula adherens maintenance.[UniProtKB/Swiss-Prot Function]



[View online »](#)