

Product datasheet for **MR210700L4V**

Zfp148 (NM_011749) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Zfp148 (NM_011749) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zfp148
Synonyms:	2210405J08Rik; AI480666; AW045217; BERF-1; BFCOL1; ZBP-89; Znf148
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011749
ORF Size:	2382 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210700).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011749.4 , NP_035879.1
RefSeq Size:	9431 bp
RefSeq ORF:	2385 bp
Locus ID:	22661
UniProt ID:	Q61624
Cytogenetics:	16 23.93 cM
Gene Summary:	Involved in transcriptional regulation. Represses the transcription of a number of genes including gastrin, stromelysin and enolase. Binds to the G-rich box in the enhancer region of these genes.[UniProtKB/Swiss-Prot Function]


[View online »](#)