

## Product datasheet for **MR210716L4V**

### Cpt1c (NM\_153679) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cpt1c (NM_153679) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cpt1c
Synonyms:	6530437J22Rik; 9630004I06Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_153679
ORF Size:	2394 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210716).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_153679.2</a> , <a href="#">NP_710146.1</a>
RefSeq Size:	2811 bp
RefSeq ORF:	2397 bp
Locus ID:	78070
UniProt ID:	<a href="#">Q8BGD5</a>
Cytogenetics:	7 B3
Gene Summary:	May play a role in lipid metabolic process.[UniProtKB/Swiss-Prot Function]



[View online »](#)