

Product datasheet for MR210730L3

Map3k20 (NM_023057) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Map3k20 (NM_023057) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Map3k20
Synonyms:	AV006891; B230120H23Rik; MLTK; MLTKalpha; MLTKbeta; Zak
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210730).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_023057
ORF Size:	2409 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_023057.5 , NP_075544.1
RefSeq Size:	3365 bp
RefSeq ORF:	2409 bp
Locus ID:	65964
UniProt ID:	Q9ESL4
Cytogenetics:	2 C3
Gene Summary:	Stress-activated component of a protein kinase signal transduction cascade. Regulates the JNK and p38 pathways (PubMed:11042189). Part of a signaling cascade that begins with the activation of the adrenergic receptor ADRA1B and leads to the activation of MAPK14 (By similarity). Pro-apoptotic. Role in regulation of S and G2 cell cycle checkpoint by direct phosphorylation of CHEK2 (By similarity). Involved in limb development (PubMed:26755636). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR210730L3