

## Product datasheet for **MR210742L4V**

### Clcn7 (NM\_011930) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Clcn7 (NM_011930) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Clcn7
Synonyms:	CIC-7; D17Wsu51e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011930
ORF Size:	2412 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210742).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011930.2</a> , <a href="#">NP_036060.1</a>
RefSeq Size:	4063 bp
RefSeq ORF:	2412 bp
Locus ID:	26373
UniProt ID:	<a href="#">O70496</a>
Cytogenetics:	17 12.53 cM
Gene Summary:	Slowly voltage-gated channel mediating the exchange of chloride ions against protons. Functions as antiporter and contributes to the acidification of the lysosome lumen. [UniProtKB/Swiss-Prot Function]



[View online »](#)