

Product datasheet for **MR210803L3V**

Nell2 (NM_016743) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nell2 (NM_016743) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nell2
Synonyms:	A330108N19Rik; mel91; R75516
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016743
ORF Size:	2460 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210803).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016743.2 , NP_058023.2
RefSeq Size:	3385 bp
RefSeq ORF:	2460 bp
Locus ID:	54003
Cytogenetics:	15 E3
Gene Summary:	Required for neuron survival through the modulation of MAPK pathways (By similarity). Involved in the regulation of hypothalamic GNRH secretion and the control of puberty (By similarity).[UniProtKB/Swiss-Prot Function]


[View online »](#)