

Product datasheet for **MR210881L4V**

Kif18b (NM_197959) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kif18b (NM_197959) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kif18b
Synonyms:	3000004C01Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_197959
ORF Size:	2502 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210881).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_197959.2 , NP_932063.2
RefSeq Size:	3345 bp
RefSeq ORF:	2505 bp
Locus ID:	70218
UniProt ID:	Q6PFD6
Cytogenetics:	11 E1
Gene Summary:	In complex with KIF2C, constitutes the major microtubule plus-end depolymerizing activity in mitotic cells. Its major role may be to transport KIF2C and/or MAPRE1 along microtubules (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)