

Product datasheet for **MR210934**

Matr3 (NM_010771) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Matr3 (NM_010771) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Matr3
Synonyms:	1110061A14Rik; 2810017I02Rik; AI841759; AW555618; D030046F20Rik; mKIAA0723
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR210934 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTCCAAGTCATTCCAGCAGTCATCTCTCGTAGGGATTCACAGGGTCATGGGCGTGACCTGTCTGCAG
 CAGGAATAGGCCTTCTTGTCTGCTACCCAGTCTTTAAGTATGCCAGCATCTCTTGAAGGATGAACCA
 GGGTACTGCACGCCTTGCTAGTTAATGAATCTTGAATGAGTCTTCATTGAATCAACAAGGAGCTCAT
 AGTGCAGTGTCTTCTGCTAGTACTTCTCCATAATTTGCAGTCTATATTTAACATTGGAAGTAGAGGTC
 CACTCCCTTTGTCTTCTCAACACCGTGGAGATACAGACCAGGCCAGTAATATTTGGCCAGCTTTGGTCT
 GTCTGCTAGAGACTTAGATGAAGTGAAGTGCATCCAGAGGACAAGATTACTCCTGAGAAGTTCGCCCAA
 ATCTTCTACAGCTTAAAAGGAGGAGAAGTGAAGAAGGCCCTACATTGAGTTATGGTAGAGATGGCAGAT
 CTGCTACACGGGAGCCACCATACAGAGTACCTAGGGATGATTGGGAAGAAAAAGGCACCTTGAAGAGA
 TAGTTTTGATGATCGTGGTCTAGTCTCAACCCAGTGTGATTATGACCATGGAAGTCGTTCTCAAGAA
 TCTGGTTATTATGACAGAATGGATTATGAAGATGACAGATTAAGAGATGGAGAAAGGTGTAGGGATGATT
 CTTTTTTGGTGAGACCTCGCATAACTATCATAAATTTGACAGTGAAGTATGAGAGAAAGGACGTGGTCC
 TGGCCCTTACAAGAGAGATCTCTTTTGAAGAAAGAGGGCGCTCCTCCAAGTAGCAATATTGAAGAC
 TTCCATGGACTCTTACCGAAGGGTTATCCCATCTGTGCTCTATATGTGATTTGCCAGTTCATTCTAATA
 AGGAGTGGAGTCAACATATCAATGGAGCAAGTCAAGTGCATGATGCCAGCTTCTTCTTGAATCTACCC
 AGAATGGAATCCTGACAATGATACTGGACACACAATGGGTGATCCATTCTGCTGCAGCAGTCTACAAAC
 CCAGCACCAGGGATTCTGGGACCGCCCTCCTTCATTTTCATCTTGGAGGACCAGCAGTTGGACCAAGAG
 GAAATCTGGGTGCTGGTAAATGGAACCTACAGGGACCAAGACACATGCAAAAAGGCAGATGGAAACCCAG
 CCGAGTTGTTTACATTATGGATTTTCAGCGAGGAAAAAACTTGAAGATACCAGCTGTACAAGTGGTAGAA
 CCATTTGGAGTCAATTTCAATCATCTGATTTTAAATAAAATTAATGAGGCATTTATTGAAATGGCCACCA
 CAGAAGATGCTCAGGCTGCTGTGGATTATTATACAACACACCAGCATTAGTGTGGCAAGCCAGTGGAG
 AGTTTCAATTTCCAGAAGTATAAGAGAATAAAGAAACCTGAAGGAAAACCCAGATCAGAAGTTTGTATCAA
 AAGCAAGAAGTGGTCTGTGATACATCTCAGCAATTTACCTCATTCTGGCTATTCTGACAGTCCCGTCC
 TCAAGCTTGTGAGCCTTATGGGAAAATAAAGAATTATATCTTAATGAGGATGAAAAGTCCAGGCTTTTAT
 TGAGATGGAACCCAGAGAAGATGCAATGGCAATGGTTGACCACTGTTAAAGAAGGCACCTTGGTTTCAG
 GGGAGATGTGTTAAAGTTGATCTGTCTGAGAAATACAAAAAACTGGTGTGCGGATTCCTAACAGAGGCA
 TTGACTTACTGAAAAAGATAAATCCCGGAAAAGATCCTATTCTCCAGATGGAAAAGAATCTCCAAGTGA
 TAAGAAGTCCAAAAGTGTGCTCAGAAGACTGAGAGTCCAGCTGAGGGTAAAGAACAAGAAGAGAAGTCCG
 GGTGAGGATGGTGAAGGACACCAAGGATGACCAGACAGAACAGGAGCCAGCATGCTTCTTGAATCTG
 AAGATGAACTGCTTGTGATGAGGAAGAAGTGCAGCACTACTAGAAAGTGGTAGTTTCAAGTGGGCGATGA
 GACCGATCTTGTAAATTTAGGGGATGTATCTTCTGATGGGAAAAAGAGCCTTCAGATAAAGCTGTGAAA
 AAAGATCCAAGTCTTCAAGCAACATCCAAGAAAAAACTTAAAAAGTGGACAAGATTGAGGAACTTGATC
 AAGAAAATGAAGCTGATTGGAAAATGGGATTAATAATGAGGAAAATACAGAACCGGTGCCGAATCTGC
 TGAGAATGCTGATGATCCAACAAGACACAAGTAAAAATGCAGATGGCCAAAATGATGAGAATAAGGAG
 GATTATACAATCCCGATGAGTATAGAATCGGACCATATCAGCCTAATGTTCTGTTGGTATAGACTATG
 TGATACCTAAAACAGGGTTTTACTGTAAGCTGTGTTCACTCTTTTATACAATGAAGAAGTTGCAAGAA
 TACTCATTGCAGCAGCCTTCTCATTATCAGAAATTAAGAAATTTCTGAATAAATTGGCAGAAGAACGG
 AGGCAGAAGAAGGAAACT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

ACCN:	NM_010771
ORF Size:	2541 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_010771.2
RefSeq Size:	4472 bp
RefSeq ORF:	2541 bp
Locus ID:	17184
UniProt ID:	Q8K310
Cytogenetics:	18 19.14 cM
MW:	94.6 kDa
Gene Summary:	May play a role in transcription or may interact with other nuclear matrix proteins to form the internal fibrogranular network. In association with the SFPQ-NONO heteromer may play a role in nuclear retention of defective RNAs (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR210934