

Product datasheet for **MR210934L4V**

Matr3 (NM_010771) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Matr3 (NM_010771) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Matr3
Synonyms:	1110061A14Rik; 2810017I02Rik; AI841759; AW555618; D030046F20Rik; mKIAA0723
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010771
ORF Size:	2541 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210934).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010771.3
RefSeq Size:	4472 bp
RefSeq ORF:	2541 bp
Locus ID:	17184
UniProt ID:	Q8K310
Cytogenetics:	18 19.14 cM
Gene Summary:	May play a role in transcription or may interact with other nuclear matrix proteins to form the internal fibrogranular network. In association with the SFPQ-NONO heteromer may play a role in nuclear retention of defective RNAs (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)