

Product datasheet for **MR210943L3V**

Arhgef15 (NM_177566) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arhgef15 (NM_177566) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arhgef15
Synonyms:	D130071N09; D530030K12Rik; E5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_177566
ORF Size:	2550 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210943).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177566.3
RefSeq Size:	4213 bp
RefSeq ORF:	2550 bp
Locus ID:	442801
UniProt ID:	Q5FWH6
Cytogenetics:	11 B3
Gene Summary:	Specific GEF for RhoA activation. Does not activate RAC1 or CDC42. Regulates vascular smooth muscle contractility. Negatively regulates excitatory synapse development by suppressing the synapse-promoting activity of EPHB2.[UniProtKB/Swiss-Prot Function]



[View online »](#)