

Product datasheet for **MR210945L3V**

Phc2 (NM_001195130) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Phc2 (NM_001195130) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Phc2
Synonyms:	A3galt2; AA415044; D4Erttd810e; D130050K19Rik; Edr2; Mph2; p36
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001195130
ORF Size:	2550 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210945).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001195130.1 , NP_001182059.1
RefSeq Size:	3949 bp
RefSeq ORF:	2553 bp
Locus ID:	54383
UniProt ID:	Q9QWH1
Cytogenetics:	4 62.12 cM



[View online »](#)

Gene Summary:

Component of a Polycomb group (PcG) multiprotein PRC1-like complex, a complex class required to maintain the transcriptionally repressive state of many genes, including Hox genes, throughout development. PcG PRC1 complex acts via chromatin remodeling and modification of histones; it mediates monoubiquitination of histone H2A 'Lys-119', rendering chromatin heritably changed in its expressibility.[UniProtKB/Swiss-Prot Function]