

Product datasheet for **MR210951L4V**

Tgfbr3 (NM_011578) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tgfbr3 (NM_011578) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tgfbr3
Synonyms:	1110036H20Rik; AU015626; AW215636; TBR111
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011578
ORF Size:	2550 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210951).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011578.3
RefSeq Size:	6074 bp
RefSeq ORF:	2553 bp
Locus ID:	21814
UniProt ID:	O88393
Cytogenetics:	5 E5
Gene Summary:	Binds to TGF-beta. Could be involved in capturing and retaining TGF-beta for presentation to the signaling receptors (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)