

Product datasheet for **MR210960L3V**

Xab2 (NM_026156) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Xab2 (NM_026156) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Xab2
Synonyms:	0610041O14Rik; AV025587
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026156
ORF Size:	2568 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210960).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026156.1 , NP_080432.1
RefSeq Size:	2654 bp
RefSeq ORF:	2568 bp
Locus ID:	67439
UniProt ID:	Q9DCD2
Cytogenetics:	8 1.92 cM
Gene Summary:	Involved in pre-mRNA splicing as component of the spliceosome. Involved in transcription-coupled repair (TCR), transcription and pre-mRNA splicing.[UniProtKB/Swiss-Prot Function]



[View online »](#)