

Product datasheet for **MR210971L4V**

Usp5 (NM_013700) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Usp5 (NM_013700) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Usp5
Synonyms:	AA407472; ISOT; ISOT-1; Ucht
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013700
ORF Size:	2574 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210971).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013700.1 , NP_038728.1
RefSeq Size:	3176 bp
RefSeq ORF:	2577 bp
Locus ID:	22225
UniProt ID:	P56399
Cytogenetics:	6 59.17 cM



[View online »](#)

Gene Summary:

Cleaves linear and branched multiubiquitin polymers with a marked preference for branched polymers. Involved in unanchored 'Lys-48'-linked polyubiquitin disassembly. Binds linear and 'Lys-63'-linked polyubiquitin with a lower affinity (By similarity).[UniProtKB/Swiss-Prot Function]