

Product datasheet for **MR210978L3V**

Cdhr1 (NM_130878) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cdhr1 (NM_130878) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cdhr1
Synonyms:	mKIAA1775; Pcdh21; Prcad
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_130878
ORF Size:	2580 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210978).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_130878.2
RefSeq Size:	4364 bp
RefSeq ORF:	2580 bp
Locus ID:	170677
UniProt ID:	Q8VHP6
Cytogenetics:	14 B
Gene Summary:	Potential calcium-dependent cell-adhesion protein. May be required for the structural integrity of the outer segment (OS) of photoreceptor cells.[UniProtKB/Swiss-Prot Function]



[View online »](#)