

Product datasheet for **MR210990**

Dpp9 (NM_172624) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dpp9 (NM_172624) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Dpp9
Synonyms:	6430584G11Rik; A330078I11; DPP IX; DPRP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Protein Sequence:

>MR210990 protein sequence
 Red=Cloning site Green=Tags(s)


MCSGVSPVEQVAAGDMDTAARFCVQKHSWDGLRSIIHGSRKSSGLIVSKAPHDFQFVQKPDESGPHSHR
 LYYLGMPYGSRENLLYSEIPKKVRKEALLLSWKQMLDHFQATPHHGVYSREEELLRERKRLGVFGITS
 YDFHSEGLFLFQASNSLFHCRDGGKNGFMVSPMKPLEIKTQCSGPRMDPKICPADPAFFSFINNSDLWV
 ANIETGEERRLTFCHQGSAGVLDNPKSAGVATFVIQEEFDRFTGCWWCPTASWEGSEGLKTLRILYEEVD
 ESEVEVIHVPSPALEERKTDYSRYPRGTGSKNPKIALKLAELQTDHQGKIVSSCEKELVQPFSSLFPKVEY
 IARAGWTRDGKYAWAMFLDRPQQRLQLVLLPPLAFIPAVESEAQRQAAARAVPKNVQPFVIYEEVTNVWI
 NVHDIHFHPPQAEGQQDFCLRANECKTGFCHLYRVTVELKTKDYDWTEPLSPTEDEFKCPIKEEVALTS
 GEWEVLSRHGSKIWNQTKLVYFQGTKDTLEHHLVVSYESAGEIVRLTTLGFSHSCSMSQSFDMFVS
 HYSSVSTPPCVHVKLSGPDDDLHKQPRFWASMMEAANCPDYVPEIFHFHTRADVQLYGMIIKPHTL
 QPGRKHPTVLFVYGGPQVQLVNNSFKGIKYLRLNTLASLGYAVVVIDGRGSCQRGLHFEGALKNQMGQVE
 IEDQVEGLQYVAEKYGFIDL SRVAIHGWSYGGFLSLMGLIHKPQVFKVAIAGAPVTWVMAYDTGYTERYM
 DYPENNQGYEAGSVALHVEKLPNEPNRLLILHGFLDENVHFHTNFLVSQLIRACKPYQLQIYPNERHS
 IRCRESGEHYEVTLLHFLQEHL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:


SgfI-MluI

Cloning Scheme:



ACCN:	NM_172624
ORF Size:	2589 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_172624.3
RefSeq Size:	3509 bp
RefSeq ORF:	2589 bp
Locus ID:	224897
UniProt ID:	Q8BVG4
Cytogenetics:	17 D
MW:	98 kDa
Gene Summary:	Dipeptidyl peptidase that cleaves off N-terminal dipeptides from proteins having a Pro or Ala residue at position 2.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR210990