

## Product datasheet for **MR210993L4V**

### **Mcm4 (NM\_008565) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Mcm4 (NM_008565) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mcm4
Synonyms:	19G; AI325074; AU045576; Cdc21; mcdc21; Mcmd4; mKIAA4003
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008565
ORF Size:	2589 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210993).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008565.2</a> , <a href="#">NP_032591.2</a>
RefSeq Size:	3589 bp
RefSeq ORF:	2589 bp
Locus ID:	17217
UniProt ID:	<a href="#">P49717</a>
Cytogenetics:	16 10.09 cM



[View online »](#)

**Gene Summary:**

Acts as component of the MCM2-7 complex (MCM complex) which is the putative replicative helicase essential for 'once per cell cycle' DNA replication initiation and elongation in eukaryotic cells. The active ATPase sites in the MCM2-7 ring are formed through the interaction surfaces of two neighboring subunits such that a critical structure of a conserved arginine finger motif is provided in trans relative to the ATP-binding site of the Walker A box of the adjacent subunit. The six ATPase active sites, however, are likely to contribute differentially to the complex helicase activity.[UniProtKB/Swiss-Prot Function]