

Product datasheet for **MR210996**

Rps6ka5 (NM_153587) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rps6ka5 (NM_153587) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rps6ka5
Synonyms:	3110005L17Rik; 6330404E13Rik; AI854034; MSK1; MSPK1; RLPK; RLSK; S6K-alpha-5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR210996 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAGGGGGAGGGCGGCCGAGCGGGCGCGGGCACCAGCGGGGACAGCGGCAGCGCGGCGAGCAGC
 TCCTTACTGTCAAGCACGAGCTGCGGACCGCCAACTTAACAGGACACGCAGAGAAGGTGGGATTGAAAA
 CTTTGAGCTGCTGAAGGTCTTAGGAAGTGGAGCTTATGGGAAAGTATTTCTAGTTCGTAATAAAGCGGC
 CATGATGCTGAAAGCTGTATGCTATGAAGTTTTAAAAAGGCAACAATAGTTCAGAAGGCTAAAACTA
 CAGAGCATACAAGGACAGAACGACAGTTCTGGAGCACATTAGACAGTCGCCATTTTGGTGACATTGCA
 CTATGCTTTCCAGACAGAGACCAAGCTTCATCTCATCTTAGATTATAAATGGTGGAGAAGCTTTTACT
 CATCTTTCTCAAAGAGAGCGTTTACAGAGCACGAAGTGCAGATCTATGTTGGAGAGATTGTGCTTGCCC
 TGGAACACCTGCACAAGTTGGGATTATATACCGTGACATTAAGCTTGAAAAATTTCTCCTTGATTCTAA
 CGCCCATGTGGTCTGACAGATTTCCGCTGAGTAAGGAGTTTGTGGCCGATGAAACTGAAAGAGCATAT
 TCTTTTGTGGAAGTATTGAATACATGGCTCCAGATATTGTCAGAGGGGGTATTTCGGGACACGACAAGG
 CAGTTGACTGGTGGAGCTTAGGTGTTCTAATGTATGAGCTACTAACTGGAGCATCTCCTTCACTGTCTGA
 TGGAGAAAAAATTTCCAGGCTGAGATATCTAGGAAATCTTAAAAAGTGAAGCTCCGATCCCCAGGAA
 ATGAGCACAGTAGCTAAAGACCTACTTCAGCGTCTCTTGTGAAAGATCCCAAGAAGCGATTGGGATGTG
 GCCACGGGATGCAGAAGAAATCAAAGAATCTCTTCTTTGAGAAGATAAAGTGGGATGACTTAGCTGC
 CAAAAAGTGCCTGCGCCATTCAAGCCGGTATCCGGGATGAGTTAGATGTGAGTAACTTCGCTGAAGAG
 TTCACAGAGATGGACCCACCTACTCTCCCGCTGCCTTGCCACAGAGCTCCGAGCGACTGTTTCAGGGCT
 ATTCCTTCGTTGCTCCATCTATTCTTCAAGCGTAACGACGCTGTTATAGATCCTTTCAGTCCATAT
 GGGAGTTGACCGTCTGGAGTGACAAATGTTGCCAGGAGTGAATGATGAAGGACTCACCTTCTATCAA
 CACTATGACCTAGATCTGAAGGACAAACCTTTGGGAGAAGGAAGTTTTTCAATTTGTCGCAAGTGTGAC
 ACAAAAAAATAACCAAGCATTTGCTGTCAAATAATCAGCAAAAGGATGGAAGCCAATACTCAGAAAGA
 AATAACAGCTCTGAAACTCTGCGAAGGACACCCCAATATTGTGAAACTACATGAAGTTTTTCAGCATCAG
 GTGGCAGCGTCTGCTCAGCCCCAGGACAGGTTGTTCTATGCTCTTCTACTCTTAGCTCTGCTTTTCA
 ACAGAAGTCTCACACGGAAGCCTGTCACCTGGACATGGCTGGTCCATAGCACCTCTCAGCTCCCTCCATT
 GCCACCTCCCATGCTGAAATTTGCTATTTATCTTCTCTGACAATGGACAGCTTCATACGTTTCTT
 TGATGGAAGTCTGAATGGCGGAGAGCTGTTGAGCGCATCAAGAGAAAGAAGCACTTCAGCGAGACTG
 AGGCCAGCTACATCATGCGGAAGCTTGTGTCAGCTGTGAGCCACATGCACGATGTTGGAGTGGTGCACAG
 AGATCTGAAGCCTGAGAATCTGTTGTTCACTGATGAAAACGACAATTTGGAAATTAAGTGAATGACTTT
 GGGTTTGACGCCTCAAGCCACCTGACAATCAGCCCTCAAGACCCCGTCTTACCCTCCACTACGCAG
 CACCGGAGCTCTTGACGCACAATGGCTACGATGAGTCTGTGACCTGTGGAGCTTGGGTGTCATCTGT
 TACAATGCTGTGAGGCGAGGTGCCATTCAGTCTCATGACAGAAGTTAACATGCACCAAGTGCAGTAGAA
 ATCATGAAGAAAATTAAGGGGAGACTTCTCCTTTGAAGGAGAAGCGTGAAGAATGTATCCCAAGAGG
 CTAAGATTTGATTCAAGGGCTTCTTACAGTTGATCCAATAAAAGGCTGAAAATGTCTGGTTTGAGATA
 CAATGAATGGCTTCAGGATGGCAGTCAGCTGCTCCTCAATCCTCTGATGACTCCAGACATCTGGGATCG
 TCAGGAGCAGCTGTGCACCTGTGTGAAAGCAACCTTTTCATGCCTTTAACAAGTATAAGAGAGAAGGAT
 TTTGCCTGCAGAATGTGATAAAGCCCCCTTTGGCTAAGAGGAGGAAAAATGAAGAGGACCAGCACCAGCAC
 AGAGACTCGCAGCAGCTCCAGCGAGAGCTCCCGCTCCTCCTCCTCCAGTCCCACGGCAAGACAACGCT
 ACCAAGACGCTGCAGCCAGCAATCCCACCGAAGGCAAGTAAACCAGACACCCTCTCCAGTTCTCAGAC

ACGCGTACGCGGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR210996 protein sequence
 Red=Cloning site Green=Tags(s)


MEGEGGSGGAGTSGDSGDGGEQLLTVKHELRTANLTGHAEKVGIENFELLKVLGTGAYGKVFVLRKISG
 HDAGKLYAMKVLKKAIVQKAKTTEHTRTERQVLEHIRQSPFLVTLHYAFQTETKLHLILDYINGGELFT
 HLSQRERFTEHEVQIYVGEIVLALHLHLGIIYRDIKLENILDSNGHVVLDFGLSKEFVADETERAY
 SFCGTIEYMAPDIVRGGDSGHDKAVDWSLGVLMYELLTGASPFVTDGEKNSQAEISRRIKLDSEPPYPQE
 MSTVAKDLLQRLLMKDPKRLGCGPRDAEEIKEHLFFEKIKWDDLAACKVPAPFKPVIRDELVDVSNFAEE
 FTEMDPTYSPAALPQSSERLFGYSFVAPSILFKRNAAVIDPLQFHMVDRPGVTNVARSSAMMKDSPFYQ
 HYDLDLKDKPLGEGSFSICRKCVMKKTNQAFVKIISKRMEANTQKEITALKLCGHPNIVKLHEVFHDQ
 VAASAQPPGVVLCSSLLALLFNRSLTRKPVTTWLVHSTSQLPPLPPMPEIVLFIILLSDNGQLHTFL
 VMELNNGGELFERIKRKHSETEASYIMRKLVSASHMHDVGVVHRDLKPENLLFDENDNLEIKVIDF
 GFARLKPPDNQPLKTPCFTLHYAAPELLTHNGYDESCDLWSLGVILYTMLSGQVPFQSHDRSLTCTSAVE
 IMKKIKKGDVFEFEAWKNVSQEAQDLIQGLLTVDPNKRLKMSGLRYNEWLQDGSQSSNPLMTPDILGS
 SGAAVHTCVKATFHAFNKYKREGFCLQNVDKAPLAKRRMKRTSTSTETRSSSSSSSSSSSQSHGKTTP
 TKTLQPSNPTEGSPDNTLQFSD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:




ACCN:	NM_153587
ORF Size:	2592 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_153587.1
RefSeq Size:	4406 bp
RefSeq ORF:	2592 bp
Locus ID:	73086
UniProt ID:	Q8C050
Cytogenetics:	12 E
MW:	96.6 kDa

Gene Summary:

Serine/threonine-protein kinase that is required for the mitogen or stress-induced phosphorylation of the transcription factors CREB1 and ATF1 and for the regulation of the transcription factors RELA, STAT3 and ETV1/ER81, and that contributes to gene activation by histone phosphorylation and functions in the regulation of inflammatory genes (By similarity) (PubMed:11553624, PubMed:11909979, PubMed:16806820). Phosphorylates CREB1 and ATF1 in response to mitogenic or stress stimuli such as UV-C irradiation, epidermal growth factor (EGF) and anisomycin (PubMed:11909979). Plays an essential role in the control of RELA transcriptional activity in response to TNF and upon glucocorticoid, associates in the cytoplasm with the glucocorticoid receptor NR3C1 and contributes to RELA inhibition and repression of inflammatory gene expression (PubMed:12628924, PubMed:16806820). In skeletal myoblasts is required for phosphorylation of RELA at 'Ser-276' during oxidative stress (PubMed:12628924). In erythropoietin-stimulated cells, is necessary for the 'Ser-727' phosphorylation of STAT3 and regulation of its transcriptional potential (PubMed:11553624). Phosphorylates ETV1/ER81 at 'Ser-191' and 'Ser-216', and thereby regulates its ability to stimulate transcription, which may be important during development and breast tumor formation (By similarity). Directly represses transcription via phosphorylation of 'Ser-1' of histone H2A (By similarity). Phosphorylates 'Ser-10' of histone H3 in response to mitogenics, stress stimuli and EGF, which results in the transcriptional activation of several immediate early genes, including proto-oncogenes c-fos/FOS and c-jun/JUN (PubMed:15870105, PubMed:16517600). May also phosphorylate 'Ser-28' of histone H3 (PubMed:11441012, PubMed:15870105). Mediates the mitogen- and stress-induced phosphorylation of high mobility group protein 1 (HMGN1/HMG14) (By similarity). In lipopolysaccharide-stimulated primary macrophages, acts downstream of the Toll-like receptor TLR4 to limit the production of pro-inflammatory cytokines (PubMed:18690222). Functions probably by inducing transcription of the MAP kinase phosphatase DUSP1 and the anti-inflammatory cytokine interleukin 10 (IL10), via CREB1 and ATF1 transcription factors (PubMed:18690222). Plays a role in neuronal cell death by mediating the downstream effects of excitotoxic injury (PubMed:12807421). Phosphorylates TRIM7 at 'Ser-106' in response to growth factor signaling via the MEK/ERK pathway, thereby stimulating its ubiquitin ligase activity (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR210996