

Product datasheet for **MR211031L4V**

OsbpI5 (BC079872) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	OsbpI5 (BC079872) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	OsbpI5
Synonyms:	ORP5, Osbp2, Obph1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC079872
ORF Size:	2622 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211031).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC079872 , AAH79872
RefSeq Size:	3834 bp
RefSeq ORF:	2624 bp
Locus ID:	79196
Cytogenetics:	7 88.29 cM



[View online »](#)

Gene Summary:

Lipid transporter involved in lipid countertransport between the endoplasmic reticulum and the plasma membrane: specifically exchanges phosphatidylserine with phosphatidylinositol 4-phosphate (PI4P), delivering phosphatidylserine to the plasma membrane in exchange for PI4P, which is degraded by the SAC1/SACM1L phosphatase in the endoplasmic reticulum. Binds phosphatidylserine and PI4P in a mutually exclusive manner. May cooperate with NPC1 to mediate the exit of cholesterol from endosomes/lysosomes. Binds 25-hydroxycholesterol and cholesterol.[UniProtKB/Swiss-Prot Function]