

Product datasheet for **MR211062L3V**

Rad54b (NM_001039556) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rad54b (NM_001039556) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rad54b
Synonyms:	E130016E03Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001039556
ORF Size:	2661 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211062).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001039556.2
RefSeq Size:	2925 bp
RefSeq ORF:	2661 bp
Locus ID:	623474
UniProt ID:	Q6PFE3
Cytogenetics:	4 A1
Gene Summary:	Transcriptional repressor that down-regulates the expression of the fibrinogen gamma chain. Represses transcription of GSK3B gene promoter via its interaction with APBA1. [UniProtKB/Swiss-Prot Function]



[View online »](#)