

Product datasheet for **MR211110L1V**

Pcdha1 (BC082326) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pcdha1 (BC082326) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pcdha1
Synonyms:	CNR
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	BC082326
ORF Size:	2706 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211110).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC082326 , AAH82326
RefSeq Size:	4942 bp
RefSeq ORF:	2708 bp
Locus ID:	116731
Cytogenetics:	18 B2-B3



[View online »](#)