

Product datasheet for **MR211141L3V**

Scube1 (BC080278) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Scube1 (BC080278) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Scube1
Synonyms:	7330410C13Rik; A630023E24Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC080278
ORF Size:	2721 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211141).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC080278.1
RefSeq Size:	5147 bp
RefSeq ORF:	2723 bp
Locus ID:	64706
Cytogenetics:	15 E1
Gene Summary:	Could function as an adhesive molecule and its matrix bound and soluble fragments may play a critical role in vascular biology.[UniProtKB/Swiss-Prot Function]



[View online »](#)