

Product datasheet for **MR211154**

Ptk2 (BC030180) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ptk2 (BC030180) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ptk2
Synonyms:	Fadk; FAK; FRNK; mKIAA4203; p125FAK
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR211154 representing BC030180.
Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCAGCTGCTTATCTTGACCCAACTTGAATCACACACCAAGTTCGAGTACTAAGACTCACCTGGGT
ACTGGCATGGAACGGTCCCCTGGTACAATGGAACGAGATTAAAGTCTTTCATTATTTGAAAGCAGT
AGTGAGCCAACCACTGGGCCAGTATTATCAGGCATGGAGATGCTACTGATGTCAGGGGCATCATTAG
AAGATAGTGGACAGTCAAAAAGTAAAGCACGTGGCCTGCTATGGATTTCCCTCAGTCACCTGCGATCA
GAGGAGGTGCACTGGCTGCACGTGGACATGGGTGTCTCCAGTGTGAGGGAGAAATATGAACGCCCCAC
CCACCAGAGGAGTGAAAATATGAGTTGAGAATTCGTTACTTGCCAAAAGGATTTCTAAACCAGTTTACT
GAAGATAAGCCAACATTGAATTTCTTCTATCAACAGGTGAAGAGTACTACATGCTAGAAATAGCTGAT
CAAGTAGACCAAGAAATAGCTTTGAAGTTGGGTTGTTGGAAATAGGCGATCCTATTGGGAGATGAGG
GGTAAATGCCCTAGAGAAGAAGTCCAACATGAAGTATTAGAAAAGATGTTGGTTTAAAGCGATTTTTT
CCTAAGAGTTTACTGGATTCTGTCAAGGCCAAAACACTAAGAAAGCTGATCCAGCAAACATTCAGGCAG
TTTGCCAACTTAATAGAGAAGAAAGCATCCTGAAATCTTTGAGATTCTTCTCCTGTGTACAGATTT
GACAAAGAGTGCTTCAAGTGTGCCCTTGGGTCAAGTTGGATCATTCTGTGGAATTGGCAATCGGCCCA
GAAGAAGGGATCAGTTACCTGACAGACAAAGGCTGCAATCCCACACACCTAGCAGACTTTAACCAAGTG
CAGACCATCCAGTACTCAAACAGTGAAGACAAAGACAGGAAAGGAATGCTACAACATAAAATCGCTGGT
GCTCCTGAGCCTTGACCGTGACGGCACCATCCCTAACCATTGCCGAGAACATGGCCGACCTGATTGAT
GGATACTGCCGGCTGGTGAATGGAGCCACACAATCTTTCATCATCAGACCCAGAAAGAAGGTGAACGG
GCGTTGCCATCAATACCAAAGTTGGCCAAAGTGAAGCAAGGCATGCCGACACATCAGTCTCTGTG
TCAGAGACAGATGACTATGCAGAGATCATCGATGAGGAAGACACATACCCATGCCCTCGACCAGGAT
TATGAGATTTCAGAGAGAAAGAATAAGAACTTGGACGCTGATTGGAGAAGGTCAGTTTGGAGATGTACAT
CAAGGCGTGTACCTGAGCCAGAGAATCCAGCTTTGGCTGTTGCAATCAAAACATGTAATAAACTGTACT
TCGGACAGCGTGAGAGAGAAGTTCCTTCAAGAAGCCTTAAACAATGCGTCAGTTTGACCATCCTCACATC
GTGAAGCTGATTGGAGTCATTACAGAGAACCCTGTCTGGATAATCATGGAGTTGTGCACACTGGAGAG
CTGAGGTCATTTTGAAGTAAGAAAATACAGCTTGGACCTGGCATCTTTGATATTATATGCCTATCAG
CTTAGTACAGCACTCGCGTATCTGGAGAGCAAAGATTTGTTACAGGGACATTGCTGCTCGGAATGTT
CTGGTGTCTTCAAATGATTGTGTAATAATGGGAGACTTTGGATTATCTCGATATATGGAAGACGACT
TACTATAAAGCTTCAAAGGAAAATTACCTATTAATGGATGGCTCCAGAGTCAATCAATTTTCGACGT
TTTACCTCAGCTAGTGACGTGTGGATGTTTGGTGTGTATGTGGGAGACTGATGCATGGCGTGAAG
CCTTTTCAAGGAGTGAAGAACAATGATGTGATCGGTGCAATTGAAAATGGGGAAAGATTACCAATGCCT
CCAAATTTGCTCTCCACCTCTACAGCCTTATGACGAAATGTTGGGCCTATGACCCAGCAGGCGGCC
AGGTTTACTGAACTAAAAGCTCAGCTCAGCACAATCCTGGAGGAGGAGAAGGTGCAGCAAGAAGAACGG
ATGAGGATGGAATCCAGAAGACAGGCTACCGTGTCTGGACTCTGGAGGGTCTGATGAAGCACCACCC
AAGCCAGCAGACCTGGTTATCCTAGCCGAGATCCAGTGAAGGATTTTATCCAGCCACAGCACATG
GTCCAGACCAATCACTACCAGTCTCTGGCTACCCGGTCCCATGGAATCCAGCCATGGCTGGCAGC
ATTTACCAGGTCAGGCATCTCTTGGACCAGACAGAATTGTGGAATCACAGACCTCAGGAGATGTCT
ATGTGGCAGCCAGTGTGGAGGATTCGGCAGCTTTGGACCTTCGAGGGATGGGACAGGTGCTTCCCCT
CACCTGATGGAGGAGCGGCTAATCCGACAGCAGGAAATGGAAGAAGACCAGCGCTGGCTGGAGAAG
GAAGAGCGATTCTGAAACCTGATGTGAGGCTCTCTCGAGGCAGCATCGACAGGGAAGACGGGAGTTTT
CAGGGTCCGACTGGAACCAACACATCTATCAGCCTGTGGGAAGCCAGATCCTGCAGCTCCACCAAAG
AAACCTCCTCGCCTGGAGCACCTGGCCACCTAAGCAACCTGTCCAGCATCAGCAGCCCTGCAGACAGC
TACAATGAGGGTGTCAAGTGGGTATCTGTGCCTGTGCAATG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
```

Protein Sequence: >Peptide sequence encoded by MR211154
 Blue=ORF Red=Cloning site Green=Tag(s)


MAAAYLDPNLNHTPSSSTKTHLGTGMERSPGTMERVLKVFHYFESSSEPTTASIIRHGDATDVRGIIQ
 KIVDSHKVKHVACYGFRLSHLRSEEVHWHVDGMGVSSVREKYELAHPPEEWKYLRIYLPKGFNLQFT
 EDKPTLNFFYQQVKS DYMLEIADQVDQEIALKLGCEIRRSYWEMRGNAL EKKS NYEVLEKDVGLKRFF
 PKSLLD SVKAKTLRKL IQQTFRQFANLNREESILKFFEILSPVYRFDKCEFCALGSSWII SVELAI GP
 EEGISYLTDKGCNPTHLADFNQVQTIQYSNSEDKDRKMLQLKIAGAPEPLTVTAPSLTIAENMADLID
 GYCVLVNGATQSFII RPQKEGERALPSIPKLANSEKQGMRTHAVSVSETDDYAEIIEDEEDTYTMPSTRD
 YEIQRERIELGRCIGEGQFGDVHQGVYLS PENPALAVAIKTCKNCTSDSVREKFLQEAL TMRQFDHPHI
 VKLIGVITENPVWIIMELCTLGELRSFLQVRKYSLDL ASLILYAYQLSTALAYLESKRFVHRDIAARNV
 LVSSNDCVKLGDFGLSRYMEDSTYYKASKGKLP IKWMAPEINFRRTSASDVWMFGVCMWEILMHGVK
 PFQGVKNDVIGRIENGERL PMPNCPPTLYSLMTKCWAYDPSRRPRFTELKAQLSTILEEEKVQQEER
 MRMESRRQATVSWDSGGSD EAPPKPSRPGYSPRSSEGFYSPQHVMVQTNHYQVSGYPGSHGIPAMAGS
 IYPGQASLLDQTELWNHRPQEMSMWQPSVEDSAALDLRGMQVLP PHLMEERLIRQQQEMEEDQRWLEK
 EERFLKPDVRLSRGSDREDDGSFQGPTGNQHIYQPVGKDP AAPPKPPRPGAPGHL SNLSSISSPADS
 YNEGKVGICACAM
 TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Recombinant protein using MR211154 also available, [TP511154](#)

Restriction Sites:

SgfI-MluI

Cloning Scheme:




ACCN:	BC030180
ORF Size:	2733 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq Size:	4681 bp
RefSeq ORF:	2735 bp
Locus ID:	14083
Cytogenetics:	15 33.94 cM
MW:	103.5 kDa

Gene Summary:

Non-receptor protein-tyrosine kinase that plays an essential role in regulating cell migration, adhesion, spreading, reorganization of the actin cytoskeleton, formation and disassembly of focal adhesions and cell protrusions, cell cycle progression, cell proliferation and apoptosis. Required for early embryonic development and placenta development. Required for embryonic angiogenesis, normal cardiomyocyte migration and proliferation, and normal heart development. Regulates axon growth and neuronal cell migration, axon branching and synapse formation; required for normal development of the nervous system. Plays a role in osteogenesis and differentiation of osteoblasts. Functions in integrin signal transduction, but also in signaling downstream of numerous growth factor receptors, G-protein coupled receptors (GPCR), EPHA2, netrin receptors and LDL receptors. Forms multisubunit signaling complexes with SRC and SRC family members upon activation; this leads to the phosphorylation of additional tyrosine residues, creating binding sites for scaffold proteins, effectors and substrates. Regulates numerous signaling pathways. Promotes activation of phosphatidylinositol 3-kinase and the AKT1 signaling cascade. Promotes activation of MAPK1/ERK2, MAPK3/ERK1 and the MAP kinase signaling cascade. Promotes localized and transient activation of guanine nucleotide exchange factors (GEFs) and GTPase-activating proteins (GAPs), and thereby modulates the activity of Rho family GTPases. Signaling via CAS family members mediates activation of RAC1. Recruits the ubiquitin ligase MDM2 to P53/TP53 in the nucleus, and thereby regulates P53/TP53 activity, P53/TP53 ubiquitination and proteasomal degradation. Phosphorylates SRC; this increases SRC kinase activity. Phosphorylates ACTN1, ARHGEF7, GRB7, RET and WASL. Promotes phosphorylation of PXN and STAT1; most likely PXN and STAT1 are phosphorylated by a SRC family kinase that is recruited to autophosphorylated PTK2/FAK1, rather than by PTK2/FAK1 itself. Promotes phosphorylation of BCAR1; GIT2 and SHC1; this requires both SRC and PTK2/FAK1. Promotes phosphorylation of BMX and PIK3R1. Isoform 9 (FRNK) does not contain a kinase domain and inhibits PTK2/FAK1 phosphorylation and signaling. Its enhanced expression can attenuate the nuclear accumulation of LPXN and limit its ability to enhance serum response factor (SRF)-dependent gene transcription (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR211154