

Product datasheet for MR211177L4V

Extl3 (NM_018788) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Extl3 (NM_018788) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Extl3
Synonyms:	2900009G18Rik; AU040105; Extl1; mKIAA0519
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018788
ORF Size:	2760 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211177).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018788.2 , NP_061258.2
RefSeq Size:	5995 bp
RefSeq ORF:	2760 bp
Locus ID:	54616
Cytogenetics:	14 D1



[View online »](#)