

## Product datasheet for **MR211228L4V**

### **Rbm10 (NM\_145627) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Rbm10 (NM_145627) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rbm10
Synonyms:	E430039K10Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145627
ORF Size:	2793 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211228).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_145627.2</a>
RefSeq Size:	3431 bp
RefSeq ORF:	2793 bp
Locus ID:	236732
UniProt ID:	<a href="#">Q99KG3</a>
Cytogenetics:	X A1.3
Gene Summary:	Not known. Binds to RNA homopolymers, with a preference for poly(G) and poly(U) and little for poly(A) (By similarity). May bind to specific miRNA hairpins (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)