

Product datasheet for **MR211279L3V**

Xrn2 (NM_011917) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Xrn2 (NM_011917) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Xrn2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011917
ORF Size:	2856 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211279).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011917.2
RefSeq Size:	3383 bp
RefSeq ORF:	2856 bp
Locus ID:	24128
UniProt ID:	Q9DBR1
Cytogenetics:	2 G2



[View online »](#)

Gene Summary:

Possesses 5'->3' exoribonuclease activity. May promote the termination of transcription by RNA polymerase II. During transcription termination, cleavage at the polyadenylation site liberates a 5' fragment which is subsequently processed to form the mature mRNA and a 3' fragment which remains attached to the elongating polymerase. The processive degradation of this 3' fragment by this protein may promote termination of transcription. Binds to RNA polymerase II (RNAP II) transcription termination R-loops formed by G-rich pause sites (By similarity).[UniProtKB/Swiss-Prot Function]